BEAUTIFLY

Glossy
PRO
series

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- www.beautifly.eu
- www.youtube.com/beautifly.eu

Device description	Input:	DC 5V
Temperature:		6-42°C (+/-5°C)
Vibration:		8000 (± 10) rpm
Charging time:		2.5 hours
Charging current:		500mA
Power:		max 8.5W
Battery Capacity:		1000mAh
Standby time:		45 days
Net weight:		145g
Device Size:		175*55*67mm

Security:

Before using the product, read this user manual and follow the recommendations contained therein when using the device.

Hot-cold principle:

Studies show that reaching the skin at a temperature of 42°C helps produce HSP protein, repair skin nerve cells and strengthen skin elasticity to prevent wrinkles. A heat treatment with a temperature of 42°C helps loosen the pores of the skin for deep cleansing, while a cold treatment with a cold effect of up to 6°C immediately tightens the pores and locks in nourishment. Alternating heat and cold care, stimulates the skin, improves circulation and accelerates the metabolism of retained fluids. In addition, it effectively eliminates dark circles under the eyes, smoothes fine lines and swelling, and brightens the skin. Combined with microvibrations, it has a positive effect on skin firmness and reduces fatigue. The skin becomes radiant and illuminated in a natural way.

EMS electrostimulation:

EMS contributes to the activation of skin cells, stimulating them to produce collagen, which improves skin elasticity, smoothes wrinkles and leaves the skin smooth and healthy.

Crio effect:

1. Promotes skin firmness, maintains elasticity, prevents allergies and other skin problems.
2. Reduces pores, retains moisture in the skin and smoothes wrinkles.
3. Regenerates muscles, relieves skin fatigue and tension, reduces redness and pain, and inhibits the multiplication of acne lesions.

Photon light therapy:

1. Red light: wavelength 620nm (±10nm) 18 LEDs deeply reach the skin 630nm to promote blood circulation, and collagen regeneration improves skin elasticity to keep the skin smooth and moisturized. Suitable for sagging and wrinkled skin.
2. Blue light: wavelength 470 (±10nm) 18 LEDs deeply reach the skin 480nm to kill bacteria, promote blood circulation, can penetrate up to 0.25mm under the skin, kill bacteria. It has anti-inflammatory and bactericidal properties, supports the treatment of acne and reduces acne scars.
3. Yellow photon light: wavelength 600nm (±10nm) 18 LEDs deeply reach the skin 610nm to improve cell oxygenation. Stimulates lymphatic drainage. Improves the condition of rough skin, used to treat skin redness.

Detailed user manual [page 2- f1]:

1. Please fully charge the device before first use. Charging time is approximately 2.5 hours. To turn on the device, press the "5" button for about 3 seconds. By default, the cleaning mode is set to intensity level 1, the device works using vibration and heat. Depending on the individual needs of the skin, the temperature can be adjusted using the "7" button. You can select photon lights by pressing the "6" button in each mode (5 minutes time for mode 1).
2. Press the Mode button to enter Cool mode (Cryotherapy), the device works with vibration and heat. According to individual needs, you can press 7 to adjust the temperature (session lasts 2 minutes. After turning off, the fan continues to run for 2 minutes to cool down)
3. Press the power button "5" to enter EMS mode. The device vibrates and produces EMS microcurrents. Depending on your individual needs, you can adjust the EMS level by

pressing the "7" button (default time setting 5 minutes for mode 3).

4. Press the "5" button to turn on the humidification mode. The device vibrates, generating heat and EMS microcurrents penetrating the skin. Depending on your individual needs, you can adjust the EMS intensity by pressing the "7" key (default time setting 5 minutes for mode 4).

5. Press the power button "5" to turn on collagen synthesis mode. The device vibrates and cools, emitting EMS microcurrent to the skin. According to your individual needs, you can adjust the EMS intensity by pressing the "7" button (default time setting 2 minutes for mode 5).

Purification Mode:

(opens pores and exports dirt from the deep layers of the skin)

Before using Glossy Pro, you should thoroughly wash off your makeup and cleanse your skin of any impurities. You can use a sonic toothbrush for this purpose. To use the Clean mode, you will need a cotton swab. Remove the transparent rim and place a cotton pad soaked in micellar fluid or other make-up removal liquid in the center of the head. Then put a transparent ring on it to attach the cotton pad. Start your device by turning on clean mode. Your skin will be heated to open the pores and better transport the dirt out, leaving it on the cotton pad. Your pores will remain clean and your skin will be ready for the next care steps.

Cryotherapy mode

(reduces and houses pores, tightens the skin)

Cool mode, i.e. Cryotherapy, effectively closes pores and maintains skin hydration. The alternating temperature from 6-42°C causes the skin to become firmer in a natural way, becoming smoother and more compact.

EMS tensioning mode

(stimulates muscle cells and tightens the skin, eliminating wrinkles)

Before starting the EMS mode, apply a serum or ampoule adapted to the needs of your skin to your face. EMS microcurrents stimulate the muscles, elasticize and energize the skin. They remove crow's feet and wrinkles on the skin surface and restore its youthful appearance, firmness and smoothness.

Humidification mode

(facilitates the absorption of nutrients into the skin, lifts and tightens)

After cleansing, apply a skin care product such as a serum, ampoule or essence to your face. Activate the moisturizing mode to transport active ingredients deep into the skin. In each mode, you can activate red, blue or yellow photon light depending on your skin's needs.

Collagen synthesis mode

(retain nourishment in the skin, activate collagen synthesis)

After absorbing skin care products, turn on the collagen synthesis mode, cool down to shrink pores and retain moisture and care in the skin. EMS microcurrents will stimulate skin cells to synthesize collagen and tighten the muscles, giving the effect of tightening and lifting the skin. In each mode, you can choose red, blue or yellow photon light depending on your skin needs.

Safety Instructions:

1. Read all safety instructions before use.
2. Follow the user manual to use this device.
3. Do not disassemble or modify the device, consult your distributor for repair.
4. Use the device only for its intended purpose, as described in this manual.
5. Do not store or use the device in a place with high temperature and humidity, such as in the sun, near fire, in the bathroom, etc. To avoid fire or electric shock.
6. Do not use or place the device in a damp condition, in case of short circuit, electric shock or fire.
7. Keep away from children and physically disabled people. Close supervision is required when the device is operated by a user under 18 years of age.
8. When the cooling function is working, the power consumption is relatively large, the fan

will have a delay in heat dissipation. The device should be set aside for approximately 10 minutes until the next use after using the Cool function.

9. Do not leave the device connected to the power supply for a long time.

10. Do not touch it with wet hands, it may result in electric shock or short circuit.

11. Do not place the device in places with high humidity, as it may result in electric shock, short circuit or insulation damage.

12. Do not leave the device within the reach of children.

13. Children must not use the device.

14. It is prohibited to use it in places where flammable materials are stored.

15. Do not clean the device with corrosive substances.

16. Use only the included power supply.

17. Before connecting the power supply, make sure the device is dry.

18. Do not use the device on sunburned skin.

19. Discontinue use in case of unwanted side effects.

Attention:

1. Stop using when the product becomes excessively hot to avoid burns.

2. Do not drop or hit the device.

3. Do not use the device if the power cord is damaged.

Contraindications:

- pregnancy and lactation
- metal implants, pacemaker
- metal objects (e.g. jewelry)
- tumor
- fever
- disease states, weakness
- disruption of the continuity of the epidermis
- previously used fillers, botox
- viral, bacterial and fungal skin diseases
- acute forms of rosacea and phlegmon
- active herpes
- epilepsy
- migraines
- thyroid diseases

Indications:

- loss of skin firmness
- gray and earthy skin color
- wrinkles
- discoloration
- dry, dehydrated skin

More information about the company and products [page 2 - f2]:

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[EN] Rules for handling electrical and electronic equipment: Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled. The device was placed on the market after August 13, 2005

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[EN] Summary declaration of conformity: The manufacturer hereby declares under its sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances.

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